



ООО «ЛАБТЕХБИО», ИНН 7806550362, ОГРН 1187847290322

КПП 780601001, ОКПО 33571120

Почтовый адрес: 190098, С-Петербург, площадь Труда, д. 4, офис 85-2.

Тел/факс+ 7 (812) 600-47-53, Моб/Viber/Whats App +7 921 918 42 70

E-mail: commerce@labtech.bio, www.labtech.bio



Ремонт оборудования ПАО.

С 2019 года наша компания является эксклюзивным дистрибьютером испанского производителя Myr (<https://myr.com.es/>), который в том числе производит известную на рынке РФ оборудование для компании Thermo : **карусельный гистопроцессор Thermo Microm STP120** и **Станцию заливки Thermo EC350**. Наши сертифицированные и обученные производителем инженеры проведут ремонт любой сложности с оригинальными запчастями, в том числе 100% восстановление до состояния нового оборудования, с гарантией 12 месяцев.

По запросу мы предоставим ссылку на завершенные сервисные контракты по 100% восстановлению приборов STP120.

СТАНЦИЯ ДЛЯ ЗАЛИВКИ ПАРАФИНОМ EC-350



**ГИСТОПРОЦЕССОР
КАРУСЕЛЬНОГО ТИПА STP120**



В условия пандемии, когда на закупку нового оборудования средства выделяются плохо, Мы предлагаем, помимо ремонта, усовершенствование гистопроцессора: устанавливая блок вытяжки испарения и подключая к вентиляции за 195 тыс. руб., устанавливая комплект удвоения производительности со 120 до 240 кассет и проводя полную модернизацию прибора за 1,115 млн. руб. с гарантией 12 месяцев.



ООО «ЛАБТЕХБИО», ИНН 7806550362, ОГРН 1187847290322

КПП 780601001, ОКПО 33571120

Почтовый адрес: 190098, С-Петербург, площадь Труда, д. 4, офис 85-2.

Тел/факс+ 7 (812) 600-47-53, Моб/Viber/Whats App +7 921 918 42 70

E-mail: commerce@labtech.bio, www.labtech.bio

Таблица. Стандартные варианты модернизации карусельного гистопроцессора STP120

Номер по каталогу	Комплектующие карусельного тканевого процессора STP-120	Комплект модернизации	Комплекс т-удвоения пр-ти	Комплекс т-экология
STP-010	Кольцо воздуховода	10 954		10 954
STP-025	Зубчатый ремень 364XL037	6 204		
STP-030	Электродвигатель 27 мин ⁻¹	48 975		
STP-032	Электродвигатель шаговый	27 730		
STP-035	Панель передняя (дисплей)	46 226		
STP-036	Плата "MICCAR V30" (цепи управления)	174 463		
STP-037	Плата "BASCAR V15" (силовые цепи)	234 373		
STP-039	Кабель соединительный плоский	6 127		
STP-040	Батарея 12v	14 828		
STP-041	Трансформатор тороидальный 8902	52 223		
STP-049	Вентилятор угольного фильтра 115V	7 677	7 677	7 677
STP-050	Угольный фильтр (2 ед. в уп.)	8 521	8 521	8 521
STP-073	Комплект очистки через угольные фильтры и системы удаления паров реагентов			103 470
STP-077	Подшипник шариковый основного вала	7 701		
STP-089	Ёмкость пластиковая для реагента (с голубыми ручками)	98 277		
STP-099	Комплект изделий для доработки до STP120-3 (2 корзины и 3 парафиновые станции)		312 849	
Запасные части, руб (Срок поставки на склад: 45 дней)		744 281	329 047	130 622
Работы, руб		372 140	164 523	65 311
Итого, руб.		1 116 421	493 570	195 933



ООО «ЛАБТЕХБИО», ИНН 7806550362, ОГРН 1187847290322

КПП 780601001, ОКПО 33571120

Почтовый адрес: 190098, С-Петербург, площадь Труда, д. 4, офис 85-2.

Тел/факс+ 7 (812) 600-47-53, Моб/Viber/Whats App +7 921 918 42 70

E-mail: commerce@labtech.bio, www.labtech.bio

Письмо об авторизации дистрибьютора Лабтехбио в России.



Letter of Authorization

Date: 01.01.2019

To whom it may concern

We, **ESPECIALIDADES MEDICAS Myr, SL** having our place of business at C/ Lleida 17-23 – 43712 Llorenç del Penedès - Spain as manufacturer of Electromedical Instruments for Histopathology and Cytology, concretely Rotary Tissue Processor model STP 120, Tissue Embedding Center models EC 350 and EC 500, Automatic Slide Stainer model SS 30 and SS 30-H and Semi-automated Rotary Microtome model M-240, all of them **Myr** branded, hereby nominate

Labtech.bio LLC.
Konogvardeyskiy blvd., building 19, office 85-2
190098, St. Petersburg
Russia

as our exclusive distributor in the territory of Russia. Labtech.bio LLC. is authorized to bid, supply, install and operational train the above-mentioned products in any tender in the assigned territory.

Labtech.bio LLC. is authorized to contact as well as correspondence and negotiate for any tender and is able to carry out maintenance and repair obligations on the product, supported by our technical expertise, as might be prescribed by tender.

This authorization letter is valid until 31.12.2020

ESPECIALIDADES MEDICAS MYR S.L.

NIF: B-45202597

C/ Lleida, 17-23. Pol. Ind. L'Empalme
43712 LLORENÇ DEL PENEDÈS - SPAIN
Telf. +34 977 668 030 Fax +34 977 668 030

Especialidades Médicas Myr, SL
Francisco Ruiz
Managing Director

Сертификат обучения .



CERTIFICATE of ACHIEVEMENT

To whom it may concern

We, Especialidades Médicas **Myr**, SL, as manufacturer of electromedical instruments for histology and cytology, hereby certify that

Tatyana Rudakova

service engineer by the company LABTECH BIO LLC., having its place of business at St. Petersburg, Russia, passed successfully the service training course on our equipment that took place in our premises during July 2019 (8th to 11th), and have acquired sufficient knowledge in all our instruments more concretely,

Embedding center models EC 350 and EC 500
Spin Tissue Processor, STP 120 and versions STP 120-1, STP 120-2 and STP 120-3
Automatic Slide Stainer versions SS 30 and SS 30-H
Semi-Automated Rotary Microtome M-240

Tatyana Rudakova has also full capacity to train end users so that they can obtain the best possible results with our **Myr** branded instruments.

And to this aim we sign present certificate in Llorenç del Penedès, July 11th, 2019

ESPECIALIDADES MEDICAS MYR S.L.

NIF. B-42262597

C/ Lleida, 17-23, Pol. Ind. L'Empalme
43712 LLORENÇ DEL PENEDÈS - SPAIN
Telf. +34 977 668 020 Fax +34 977 668 030

Especialidades Médicas **Myr**, SL

Francisco Ruiz

Managing Director

Especialidades Médicas **Myr** S.L.

C/ Lleida 17-23, Pol. Ind. L'Empalme - 43712 Llorenç del Penedès - (SPAIN)

Phone: + 34 977 66 80 20 - Fax: +34 977 66 80 30 - esp.medicas@myr.com.es - www.myr.com.es



Myr

Especialidades Médicas Myr S.L.



CERTIFICATE of ACHIEVEMENT

To whom it may concern:

We, Especialidades Médicas Myr, SL, as manufacturer of electromedical instruments for histology and cytology, hereby certify that

Valeri Iakolev

service engineer by the company LABTECH BIO LLC., having its place of business at St. Petersburg, Russia, passed successfully the service training course on our equipment that took place in our premises during July 2019 (8th to 11th), and have acquired sufficient knowledge in all our instruments more concretely,

Embedding center models EC 350 and EC 500

Spin Tissue Processor, STP 120 and versions STP 120-1, STP 120-2 and STP 120-3

Automatic Slide Stainer versions SS 30 and SS 30-H

Semi-Automated Rotary Microtome M-240

Valeri Iakolev has also full capacity to train end users so that they can obtain the best possible results with our Myr branded instruments.

And to this aim we sign present certificate in Llorenç del Penedès, July 11th, 2019

ESPECIALIDADES MEDICAS MYR S.L.

NIF: B-22022067

C/ Lluís, 17-18 - Poble Nou, L'Empordà

43712 LLORENÇ DEL PENEDÈS - SPAIN

Tel: +34 977 668 020 - Fax +34 977 668 020

Especialidades Médicas Myr, SL

Francisco Ruiz

Managing Director